

Characterization Of Antibody-Drug Conjugates Using HDX and Limited Proteolysis-Mass Spectrometry

Liana Tsiatsiani, Yunus Sariçay, Lang Xin Zhang and Eef H.C. Dirksen


Byondis B.V., Microweg 22, 6545 CM, Nijmegen, the Netherlands

BACKGROUND


The efficacy of biopharmaceuticals, such as monoclonal antibody-drug conjugates (ADC), is strongly dependent on maintaining structural integrity throughout industrial process manufacturing, and this is traditionally monitored using **spectroscopic techniques**. Such methods monitor important aspects of protein structure but offer limited resolution regarding the exact sites of structural change. At the same time, **mass spectrometry (MS)-based methods** are gaining ground in the pharmaceutical industry because they allow peptides indicative of structural integrity to be detected with shotgun or bottom-up proteomics workflows in a single chromatographic run, and thus localize structural changes with increased accuracy. In this work, we set out to evaluate the complementarity of **limited proteolysis-MS (LiP)** and **hydrogen-deuterium exchange-MS (HDX)** for the analysis of protein **higher-order structure (HOS)** in a biopharmaceutical setup.

METHODS


HDX-MS workflow




STUDY DESIGN




LiP-MS workflow



RESULTS



HDX MS showed very precisely regions that have undergone conformational changes upon linker-drug conjugation, as compared to the non-conjugated species, DAR0. The differences in solvent exposure between the different DAR species were small, but the acquired spatiotemporal deuterium uptake allowed us to locate precisely the affected sites in **VL, CL, CH1 and CH2 domains** with near amino acid resolution.





CONCLUSION

MS-based assays **confirm minor structural changes** resulting from drug conjugation and **locate** the exact sites in CH1 and CL domains. The **choice** of structural MS-method depends on the **purpose** of the investigation.

Take-home message


Both spectroscopic and MS-based methods are **sensitive** enough to detect subtle changes in protein HOS but they **differ** greatly in **resolution** and **specificity**.

Chymotryptic LiP-MS heat map



A dual enzyme LiP-MS approach revealed regions in CL and CH1 domains that have become more accessible for proteolysis upon drug conjugation. These regions are located near the conjugation sites.

Tryptic LiP-MS heat map



x17-fold